



Same Day Shipping On Most Items!



Our Mission



Preferred Medical Products LLC

is dedicated to providing the health care industry with quality products, while maintaining affordable medical care to our greatest asset
...the **PATIENT!**

Founded in 1991, Preferred Medical Products is
an American owned and operated company offering the highest
quality sterile and non-sterile equipment covers to the medical industry.

GENERAL INFORMATION

Mail: PO Box 100
Physical: 191 Industrial Dr.
Ducktown, TN 37326

Phone: (423) 496-3738
(800) 441-1161

Fax: (423) 496-3758

Web: www.preferred-medical.com



EC	REP
----	-----

mdi EUROPA GmbH
Langenhagener Str. 71
D-30855 Langenhagen
Germany
Phone +49-511-39 08 95 30
Fax +49-511-39 08 95 39

Bag Decanters



Lead Protection



Lead Glasses • Lead Aprons • Sorbx®



thru page 20



13



Tablet & Phone Covers



11

Specialty Products



Shipping Policies

24

Ordering Information



23

Equipment Covers

Camera Covers

10

C-arm Drapes

12

Cassette Covers

6

Cover Kaps

6

Eazy Covers

5

Fluoro Drapes

8

Probe Covers

7

Set-Up Cover

8

Bag Decanters



PMP's BAG DECANTER is an economical, single use product designed to facilitate the aseptic transfer of fluids from flexible IV containers. The decanter's spout can be resealed to stop the flow of the IV bag's liquid contents as needed. The product is easy to use in four simple steps.
(See below)

Part #: BD-100
50 per case
(1 case minimum)

**Call our
Sales
Department
for
FREE
SAMPLES!
800-441-1161**



- 1.) Firmly grasp the decanter's container pouch at the end and peel open to expose the spiked end of the Bag Decanter.
- 2.) Hold the product (still in the container pouch) and remove the protective cap from the spiked end.
- 3.) Firmly insert the spike into the IV bag outlet, then remove the packaging from the spout.
- 4.) Remove the spout cap and squeeze the container to dispense the contents as needed.



Eazy Covers

Eazy Covers, with their elastic band closure, are easily applied. Also available without bands in most sizes. They are single-use equipment covers serving as effective bacterial and contagion barriers for a wide variety of medical equipment. All Eazy Covers are latex free and economical, costing from 25%-50% less than most comparable products available.

Product No.	Dimensions	Qty/Case
EZ-0507	(5" X 7")	50
EZ-0512	(5" X 12")	50
EZ-0548-DB	(5" X 48")	20
EZ-0572	(5" X 72")	20
EZ-0709	(7" X 9")	50
EZ-0912	(9" X 12")	50
EZ-1060	(10" X 60")	20
EZ-1208	(12" X 8")	20
EZ-1212	(12" X 12")	20
EZ-1515	(15" X 15")	20
EZ-1636	(16" X 36")	20
EZ-2012	(20" X 12")	20
EZ-2015	(20" X 15")	20
EZ-2020	(20" X 20")	20
EZ-2030	(20" X 30")	20
EZ-2226	(22" X 26")	20
EZ-2236	(22" X 36")	20
EZ-2458	(24" X 58")	20
EZ-2515	(25" X 15")	20
EZ-2640	(26" X 40")	20
EZ-2650	(26" X 50")	20
EZ-12	(28" X 12")	20
EZ-3015	(30" X 15")	20
EZ-3018	(30" X 18")	20
EZ-3020	(30" X 20")	20
EZ-3025	(30" X 25")	20
EZ-3030	(30" X 30")	20
EZ-3036	(30" X 36")	20
EZ-3040	(30" X 40")	20
EZ-3612	(36" X 12")	20
EZ-3616	(36" X 16")	20
EZ-20	(36" X 20")	20
EZ-3626	(36" X 26")	20
EZ-28	(36" X 28")	20
EZ-3630	(36" X 30")	20
EZ-3636	(36" X 36")	20
EZ-3640	(36" X 40")	20
EZ-3654	(36" X 54")	20

Product No.	Dimensions	Qty/Case
EZ-3690	(36" X 90")	20
EZ-4020	(40" X 20")	20
EZ-4028	(40" X 28")	20
EZ-4030	(40" X 30")	20
EZ-4040	(40" X 40")	20
EZ-4060	(40" X 60")	20
EZ-4220	(42" X 20")	20
EZ-4230	(42" X 30")	20
EZ-4236	(42" X 36")	20
EZ-4425	(44" X 25")	20
EZ-4440	(44" X 40")	20
EZ-4448	(44" X 48")	20
EZ-4525	(45" X 25")	20
EZ-4646	(46" X 46")	20
EZ-4836	(48" X 36")	20
EZ-4865	(48" X 65")	20
EZ-5040	(50" X 40")	20
EZ-6030	(60" X 30")	20
RB-4060	(40" X 60") Telescopic Fold	20

*These Eazy Covers are available with adhesive tabs for an additional \$.15 per unit

Custom Sizes Available!!!



Some covers available in space saving 100 count boxes!

All sizes of Eazy Covers are also available Non Sterile!

Ask for availability!

Cover Kaps

Product	Dimensions	Qty/Case
SC-1519	11" depth	20
SC-20	9" depth	20
SC-25	13" depth	20
SC-30	14" depth	20
SC-35	18" depth	20
SC-40	20" depth	20
SC-45	22" depth	20
SC-60	27" depth	20
SC-70	31" depth	20
SC-80	35" depth	20



Cover Kaps are economical, easily applied, single-use products that function as effective bacterial and contagion barriers for assorted medical equipment. The Kap's domed shape and elastic band closure facilitates a wide variety of applications.



Product	Dimensions	Qty/Case
SC-4216	42" x 16" (sewn corners)	20
SC-5416	54" x 16" (sewn corners)	20
SC-6020	60" x 20" (sewn corners)	20
SC-7624	76" x 24" (sewn corners)	20



Cassette Covers

Non-Sterile
Available



Product	Dimensions	Qty/Case
CC-16.523	(16.5" X 23")	20
CC-1620 W/S	(16" X 20") Removable Sleeve	20
CC-2420 W/S	(24" X 20") Removable Sleeve	20
CC-1820	(18" X 20")	20
CC-2040	(20" x 40")	20
CC-2040 W/T	(20" X 40") W/Adhesive	20
CC-2324.5	(23" X 24.5")	20
CC-2324.5 W/T	(23" X 24.5") W/Adhesive	20
CC-2426 W/T	(24" X 26") W/Adhesive	20
CC-2436 W/T	(24" X 36") W/Adhesive	20
CC-3050	(30" X 50")	20

Probe Covers



Product	Dimensions	Case	Bands	Adhesive	Gel	Straight End	Tapered End	Scrunch Fold	Accordion Fold	Telescopic Fold	Straight Fold
PC-1.3836-WS	1.38" x 36" with sheath	50									
PC-0314 DLT	3" X 14"	20		•			•				•
PC-0332 SB	3" X 32"	20	•	•		•					•
PC-0412 SB	4" X 12"	50	•			•					•
PC-0433 SB	4" X 33"	50	•			•					•
PC-0433 STG	4" X 33"	50		•	•	•					•
PC-0520 W/B	5" X 20"	20	•				•				•
PC-0572-TF	5" X 72"	20		•						•	
PC-0624 W/B	6" X 24"	20	•				•	•			
PC-0648-AF	6" X 48"	20	•						•		
PC-0648 A1G	6" X 48" (1 mil plastic)	20	•		•				•		
PC-0648 W/B	6" X 48"	20	•				•	•			
PC-0648 W/B1	6" X 48" (1 mil plastic)	20	•				•	•			
PC-0696 W/B	6" X 96"	20	•				•	•			
PC-0760-TF	7" X 60"	20		•						•	
PC-0872-TF	8" X 72"	20		•						•	
PC-0896-TF	8" X 96"	20		•						•	
PC-0896 W/B	8" X 96"	20	•				•	•			
PC-1072-TF	10" X 72"	20								•	



PC-0314 DLT used on gamma finder

**Sterile
Gel**

Available in
20 count boxes

Set-Up Covers



Pictured is product PD-5490

Set-Up Covers are useful for isolating tables, trays and other horizontal surfaces where adhesion is not necessary.

Product	Dimensions	Qty/Case
PD-2754	27" X 54"	20
PD-4470	44" X 70"	20
PD-4482	44" X 82"	20
PD-4482-NS	44" X 82"	50
PD-4545	45" X 45"	20
PD-5454	54" X 54"	20
PD-5470	54" X 70"	20
PD-5490	54" X 90"	20

Non Sterile represented by -NS at the end of product name



Fluoro Drapes

Fluoro Drapes have adhesive strips (with an easy-peel protective backing) along the edges, with which the drape may be firmly affixed. This adhesive is formulated to leave no residue and allows easy removal after use.

Product	Dimensions	Qty/Case
FD-1824	18" X 24"	20
FD-2217	22" X 17"	20
FD-2420	24" X 20"	20
FD-2430-2	24" X 30"	20
FD-3620	36" X 20"	20
FD-3620-1/S	36" X 20"	20
FD-3636	36" X 36"	20
FD-3640	36" X 40"	20
FD-3644	36" X 44"	20
FD-3660-1	36" X 60"	20
FD-4444	44" X 44"	20
FD-4460	44" X 60"	20
FD-4460-1/S	44" X 60"	20



Don't want to settle for a generic cover?

We take pride in our ability to

CUSTOMIZE

Customization is PMPs forte, we love a challenge. Preferred Medical Products has built a strong reputation in the medical field as the company that will do what is necessary to provide the best fitting cover for the customer's application. We have the capability and are willing to work with our customers to assure that they receive the type cover they desire. All products are manufactured on site which gives us and our customers many options.

≧Let us design your custom cover and provide you with FREE samples≦

STERILE or NON-STERILE



★ 100% American Made ★


Camera Covers

Camera Covers are economical, telescopically folded disposable products designed to provide an effective bacterial and contagion barrier for laser camera equipment.




Product	Dimensions	Qty/Case
CAMC-0448-TF	4" X 48" Telescopic Fold	20
CAMC-0596-P	5" X 96" Perforated	20
CAMC-0596-TF	5" X 96" Telescopic Fold	20
CAMC-0796-TF	7" X 96" Telescopic Fold	20
CAMC-1096-TF	10" X 96" Telescopic Fold	20

Sponge Counters



Build your custom sponge counter using the system pictured.




*EO includes blue wire ties that allow pockets to remain open for easier insertion of sponge (pictured)

Mayo Stand Cover

Product	Dimensions	Qty/Case
MS-2454	24" X 54"	20



Sterile Mobile Device Covers



**Don't Miss
URGENT
Information!**

Protect your electronic devices during procedures with PMPs new line of sterile covers for tablets and cell phones!



Covers available for all devices with and without cases



Stay connected
in the OR!



Optional velcro for
attaching to work area

Call for free samples!
800-441-1161

Sold in cases of 50 each * Also available non-sterile

C-arm Drapes

C-arm Draping Systems are economical, easily applied and disposable. They are designed to provide an effective bacterial and contagion barrier for C-Arm X-ray equipment.



Pictured at top left and bottom left is product C-arm 6085; right is C-arm 90B.

C-arm Draping Systems


Product	Dimensions	Qty/Case
C-ARM 54	54" Tail	20
C-ARM 4060	40" X 60" (Bag w/straps)	20
C-ARM 4072	40" X 72"	20
C-ARM 4174	41" X 74" (Bag w/straps)	20
C-ARM 4265	42" X 65" (Bag w/adhesive)	20
C-ARM 4470	44" X 70" (Bag w/straps)	20
C-ARM 6085	60" X 85" (Bag w/straps)	20
C-ARM 4065B	40" X 65" (Bag w/band & straps)	20
C-ARM 54B	36" X 54" (Banded bag w/extra band)	20
C-ARM 90B	36" X 90" (Banded bag w/extra bands)	20

Mini C-Arm Draping System (U-Shaped Style)

Product	Dimensions	Qty/Case
C-ARM 4264U	42" X 64" (Strap Tie)	20
C-ARM 4885M	48" X 85" (Band & Clips)	20
C-ARM 70M	44" X 70" (Band & Clips)	20

Are You Protected?

THE PREFERRED SOLUTION TO RADIATION PROTECTION



Preferred Medical Products offers a wide variety of radiation protection for physicians and clinicians. With Interventional procedures producing extremely harmful scatter, proper protection is a must. Sorbx® used in conjunction with standard radiation protection gives you the best line of



defense. Aprons are a great start to protecting your biggest asset “YOU”, but you can’t stop there. We have learned there are better methods of protection. You can use glasses and thyroid collars and even lead gloves with your aprons, but what about your unprotected areas? With Sorbx® you are reducing the exposure to the protected areas and also providing your only line of defense to the unprotected areas. Have you ever thought of how you’re protecting your head? On the following pages you will see additional information about Sorbx® and what it can do for you and your team.

Sorbx®

Multi-Purpose Shields

Interventionalist radiation exposure is primarily from scatter radiation from the patient

Allowable Radiation Dosage (Rads) Per Year




Available without fenestration or. . . .



These shields are designed with an overall dimension of 12" x 16.5" and do not have a fenestration punch area. These shields are made without fenestrations for greater precision when directing the radiation beam.

The Multi-Purpose Shields are designed for use in biopsy and biliary procedures. These procedures include shunt procedures, biliary drainage procedures, spinal procedures, stents and diagnostic procedures. These shields can be used anywhere on the body depending on the entry point.

. . . .with fenestration



The shields with fenestration are designed with an overall dimension of 12" x 16.5" and have a 3.5" x 2" rectangular fenestration area. The SorbX® protection area is 12" x 15.5". These shields are designed for use in various biopsy and nephrostomy (kidney) procedures. These shields can be used anywhere on the body depending on the entry point.

The SorbX® Dialysis Drape has two layers of radiation protection with an overall dimension of 38.5" x 36". It is designed with a 4.5" x 9" rectangular fenestration area for needle or catheter access, and includes a 23" x 32" absorbent area. Velcro strips containing SorbX® are included allowing the doctor to reposition them based on the angle of the image they are taking.


The Dialysis Drape is used on dialysis access and declotting procedures. The access point is the underside of the forearm and the needle or catheter is placed into the fistula to remove blood clots.



At 90 KvP	1 Layer	2 Layers
Attenuation	75.3%	90.8%
Lead Equivalency (mm Pb)	.116	.278

BIODEGRADABLE RADIATION PROTECTION

Provides clinicians with additional protection from dangerous scatter radiation which allows them to perform more procedures while increasing their safety from dangerous ionized radiation.



Customization Available

Procedure specific designs combined with proper positioning allow for the maximum shaded area of protection to be formed

More about SorbX®

- **Sterile Products** – Single use sterile product
- **Green Products** – Friendly to users of the product and safe for disposal
- **Disposable** – Product can be disposed of with regular hospital waste
- **Lead-Free** – The core material is made from a proprietary blend of fine metal powders allowing it to be 100% lead free
- **Environmentally Friendly** – Radiation material is biodegradable
- **Flexible** – the SorbX® product can be folded over and will hold in position
- **Lightweight** – SorbX® is half the weight of some competitive product

OUTER LAYER ABSORBENCY

If fluid is not absorbed within a few seconds, it will get on the patient, compromise the sterile field, and run-off the table on to the floor. When shields are placed on the patients they are never flat, but are contoured based on the shape of the patients body. This places further importance on the rate fluid is absorbed. We have tested the absorbency of SorbX® material to be up to eight times faster than the competition's.

Product #	Description	Qty/CS
SX8100	Single Layer Multi-Purpose Shield, 12" x 16.5", Absorbent, No Fenestration	15
SX8110	Single Layer Multi-Purpose Shield, 12" x 16.5", Absorbent, Rectangular Fenestration	15
SX8200	Two Layer Multi-Purpose Shield, 12" x 16.5", Absorbent, No Fenestration	15
SX8210	Two Layer Multi-Purpose Shield, 12" x 16.5", Absorbent, Rectangular Fenestration	15
SX7140	Dialysis Access Declotting Drape, 38.5" x 36" (Fenestration 5 x 9)	15



Leaded Eyewear

All Front Lenses Are .75mm

All Side Shields Are .5mm (add \$15⁰⁰)

Allow a 2 - 3 week lead time for prescription glasses!!



“Interventional radiologists are at high risk for radiation-induced eye injury.”



FALCON

Adjustable nose piece and spring hinge
Available in tortoise and silver
All scripts



RAFTER

Adjustable temple length and frame angle
All scripts available.
Splash protection on top
Shown with side shields



KAYAK

Nylon frames allow maximum flexibility in fit
Large molded side splash protection
All scripts
Shown with side shields



PREFERRED LITE

With and without side shields - No Scripts
The curved front and temples with rubber tips, assures a close secure fit.
Black, crimson, bronze, purple, silver



COHUTTA

Snug fitting, curved frame
Available in blue and black
With or without side shields
No Scripts



TOCCOA

Designed to fit over prescription eyewear
Use as a goggle for splash protection
Black or tortoise



CHEROKEE

Attractive two tone temple with curved frame
No scripts, no side shields
Smaller lens reduces weight
Black, grey, silver, blue



BRAVES 2

Rubber nose pads and temple ends
With and without side shields
Most single and bifocal scripts
Silver, blue & black



OCOEE

Ultra soft fixed nose piece for added comfort and closer fit
Most scripts possible
With and without side shields
Crimson, black, silver, blue

One year warranty
against any
manufacturing
defects

- ☒ Current prescription needed for prescriptions
- ☒ Prescription needs to be reviewed prior to final selection of frame.
- ☒ Most are available for same day shipping. Please call for availability.

All glasses shipped in protective case with cloth croakie support straps to help secure the glasses to the face

Front Protection Deluxe Aprons

SELECT THE PROTECTIVE MATERIAL AND APRON STYLE THAT FITS YOUR NEEDS

Please allow a 2-3 week lead time!

REGULAR:

This is the traditional leaded protective material. The protection is .5mm Pb equiv. at 100 KVP. This is the lowest cost but the heaviest. A good product if the user is only using it for a short period of time.

FEATHER LITE:

This is the lightest leaded product on the market today. It also provides .5mmPb equiv. at 100 KVP.

ECOLITE:

This product incorporates a 2-element attenuating material and DOES NOT include lead. EcoLite is environmentally safe and provides .5mmPb equiv. at 100 KVP. You can't get lighter.

TIE APRON

MALE		FEMALE
SMT	Small	SFT
MMT	Medium	MFT
LMT	Large	LFT
XMT	Xtra-Large	XFT

DROP AWAY

MALE		FEMALE
SMDA	Small	SFDA
MMDA	Medium	MFDA
LMDA	Large	LFDA
XMDA	Xtra-Large	XFDA

BUCKLE

MALE		FEMALE
SMB	Small	SFB
MMB	Medium	MFB
LMB	Large	LFB
XMB	Xtra-Large	XFB

VELCRO

MALE		FEMALE
SMV	Small	SFV
MMV	Medium	MFV
LMV	Large	LFV
XMV	Xtra-Large	XFV

CUMMERBUND

MALE		FEMALE
SMC	Small	SFC
MMC	Medium	MFC
LMC	Large	LFC
XMC	Xtra-Large	XFC

STRETCHY BACK

MALE		FEMALE
SMSB	Small	SFSB
MMSB	Medium	MFSB
LMSB	Large	LFBS
XMSB	Xtra-Large	XFBS

All Front Protection .5MMPB EQUIV. AT 100KVP

CALL YOUR SALES REP FOR MORE INFORMATION!!

Full Protection Deluxe Aprons

Leaded (*Feather Lite and Regular*) and Lead Free (*EcoLite*)
Vest and Jacket Front Panels .5MM/Back Panels .25MM



VEST/SKIRT

MALE		FEMALE
SMTP	Small	SFTP
MMTP	Medium	MFTP
LMTP	Large	LFTP
XMTP	Xtra-Large	XFTP



2PIECE OVERLAP

MALE		FEMALE
SMTPOL	Small	SFTPOL
MMTPOL	Medium	MFTPOL
LMTPOL	Large	LFTP OL
XMTPOL	Xtra-Large	XFTPOL



OVERLAP JACKET

MALE		FEMALE
SMJOL	Small	SFJOL
MMJOL	Medium	MFJOL
LMJOL	Large	LFJOL
XMJOL	Xtra-Large	XFJOL

Over lap aprons have .25mm front panels and overlap to provide .5mm protection. They are lighter than the regular jacket and two piece aprons.

THYROID COLLARS



Tethered to Apron
TS-5



Separate Collar
TS-1



JACKET

MALE		FEMALE
SMJ	Small	SFJ
MMJ	Medium	MFJ
LMJ	Large	LFJ
XMJ	Xtra-Large	XFJ

Gonapron™ and Covers Available



PRODUCT #:
Gonapron



PRODUCT #:
NB-2016 W/T
NB-2019 W/TP

SIZING CHART

MALE							FEMALE						
TWO PIECE													
VEST				SKIRT			VEST				SKIRT		
SIZE	HEIGHT	CHEST	LENGTH	WAIST	HIPS	LENGTH	SIZE	HEIGHT	CHEST	LENGTH	WAIST	HIPS	LENGTH
S	5'5"	39"	22"	28"	31"	20"	S	5'2"	33"	19"	28"	34"	20"
M	5'7"	42"	24"	34"	38"	20"	M	5'5"	36"	20"	31"	38"	22.5"
L	5'11"	45"	25"	39"	44"	21"	L	5'8"	43"	22"	34"	41"	23"
XLG	6'2"	50"	26"	42"	46"	21"	XLG	5'10"	46"	22"	37"	44"	23"
FRONT APRON													
SIZE	HEIGHT	CHEST	LENGTH	WIDTH	JACKET SIZE		SIZE	HEIGHT	CHEST	LENGTH	WIDTH	DRESS SIZE	
S	5'5"	38"	34"	22.5"	38		S	5'2"	34"	34"	21"	4 - 6	
M	5'7"	40-43"	37"	23"	42		M	5'5"	36-39"	36"	22"	8 - 10	
L	5'11"	44-46"	40.5"	24"	46		L	5'7"	40-44"	38.5"	24"	12 - 14	
XLG	6'2"	50-54"	42.5"	26"	50		XLG	5'11"	46-50"	40.5"	26"	16 - 18	
XXLG		54-58"			54		XXLG		51-54"	42.5"			
JACKET													
SIZE	HEIGHT	CHEST	WAIST	HIPS	LENGTH		SIZE	HEIGHT	CHEST	HIPS	LENGTH		
S	5'5"	40"	30"	40"	36"		S	5'2"	36"	38"	34"		
M	5'7"	43"	37"	43"	38"		M	5'5"	38"	42"	37"		
L	5'10"	45"	42"	46"	40.5"		L	5'8"	42"	46"	40.5"		
XLG	5'11"	50" MAX	50"	52"	42.5"		XLG	5'10"	46"	50"	42.5"		


Lead Apron Racks



Two Peg Rack

Item #5139

- Holds one apron, any style
- Two mounting holes
- High tensile aluminum



Eight Peg Rack

Item #5136

- Designed to hold four aprons
- Four mounting holes
- High tensile aluminum
- 24" long



Valet Rack

Item #9890

- H 54" W 35" D 24"
- Holds up to twelve aprons
- Almond coated metal
- Hangers ordered separately



Porta Rack

Item #9880

- Star shaped base - prevents tipping
- Totally portable
- Holds Eight Aprons

STEPS IN ORDERING

STANDARD APRON WEIGHT CHART

APRON STYLE	SIZE	REGULAR	FEATHER	ECO
FRONT APRONS	S	8 LBS	7 LBS	6 LBS
	M	8.8 LBS	7.8 LBS	6.8 LBS
	L	9.6 LBS	8.5 LBS	7.5 LBS
VEST & SKIRT	S	12 LBS	11 LBS	10 LBS
	M	14.3 LBS	13 LBS	12 LBS
	L	15 LBS	14 LBS	13 LBS
JACKET	S	11.5 LBS	10.5 LBS	9.5 LBS
	M	13.5 LBS	12.5 LBS	11.5 LBS
	L	16.8 LBS	15.8 LBS	14.8 LBS
VEST & SKIRT OVERLAP	S	10 LBS	9 LBS	8 LBS
	M	12 LBS	10.5 LBS	9.5 LBS
	L	11.5 LBS	11.5 LBS	10.5 LBS
JACKET OVERLAP	S	12 LBS	11 LBS	10 LBS
	M	10 LBS	9 LBS	8 LBS
	L	13 LBS	12 LBS	11 LBS

**WEIGHT (+ OR - .5 LB)

A PROTECTIVE APRON

- 1) Choose a STYLE apron that suits your needs.
All deluxe aprons include shoulder pads and pockets.
- 2) ECONOMY APRONS are available in some styles without shoulder pads or pockets.
- 3) Determine the PROTECTIVE MATERIAL you want**
- 4) Use sizing chart to pick SIZE
- 5) Use appropriate CATALOG NUMBER
- 6) Choose a FABRIC color. We will use black on the inside of the thyroid collar, the apron and the binding
- 7) We can EMBROIDER the pocket with 20 characters.
Each line of embroidery is \$8.00. Each logo is \$10.00.

FABRIC SWATCHES




Ripstop pattern magnified

Regular



Camo



Designer



Ripstop



Safari**



Sparkles**



Galaxy**



Diamonds**



Swirls**



**Bottom row available for additional cost

Pattern on swatches is smaller than appears. Swatch size is 2" x 4"

Additional Lead Accessories may be available upon request!

Specialty Products

If you are not able to find the style or size that you require,
keep in mind that we provide custom items.

Call your sales representative for more information!

All Items Are
Latex Free

Product	Dimensions	Qty/Case
ARMBOARD COVERS		
NB-0830-C	8" X 30" w/cuff	20
INSTRUMENT POUCHES		
IP-1011	10" X 11" Divided Pouch	50
IP-1012	10" X 12" (1 lg pouch)	50
IP-2424-L	24" X 24" w/lip	20
NEEDLE BIOPSY GUN COVER (Non-Disposable Triangular Sleeve)		
PC-1000	6" X 11"	50
AUTOMATIC PULLBACK DEVICE COVER		
PC-1072-TF	10" X 72"	20
IV POLE COVER		
EZ-1636	16" X 36"	20
TRASH BAGS		
TB-1820 W/T	18" X 20" w/tape	100
MONITOR COVERS		
MC-2215-NS	22" X 15"	50
MC-2727-NS	27" X 27"	50
MC-2754-NS	27" X 54"	50
SC-5416	54" X 16" (sewn corners)	20
SC-7624	76" X 24" (sewn corners)	20
Keyboard Covers (with adhesive)		
KC-2016-NS W/T	20" X 16"	50
KC-2418-NS W/T	24" X 18"	50
KC-3024-NS W/T	30" X 24"	50

REMOTE CONTROL COVERS (Adhesive Closure)		
Product	Dimensions	Qty/Case
RMC-0410 W/T	4" X 10"	50
RMC-0512 W/T	5" X 12"	50
RMC-0814 W/T	8" X 14"	50



For the most accurate fit call with your dimensions.
Your Sales Rep can send you custom samples!

Free
Samples
Available

Non Sterile

Non Sterile items
are available
sterile also!




Used as various equipment, foot
pedal and foot switch covers

To order call toll free
800-441-1161

or fax to
423-496-3758

Stereotactic Table Absorbent Pads - 50 p/case

Product	Dimensions
AP-1220-NS	12" X 20"
AP-3648-NS	36" X 48"



Pictured (L) AP-1220-NS (C) AP-3648-NS (R) AP-3648-NS used on floor to absorb fluid

Foot Pedal Covers

Product	Dimensions	Qty/Case
EZ-0709-NS	7" X 9"	50
EZ-0912-NS	9" X 12"	50
EZ-1212-NS	12" X 12"	50
EZ-1426-NS	14" X 26"	50
EZ-1515-NS	15" X 15"	50
EZ-2012-NS	20" X 12"	50
EZ-2015-NS	20" X 15"	50
EZ-2020-NS	20" X 20"	50
EZ-2030-NS	20" X 30"	50

CT Head Rest Covers - 50 p/case

Product	Dimensions
HR-1812-NS	18" X 12"
EZ-2015-NS	20" X 15"

Cassette Covers - 100 p/case (**50 per case)

Product	Dimensions
CC-1012-L-NS	10"X12" w/lip
CC-16.523-L-NS	16.5" X 23" w/lip
CC-1820-NS	18" X 20"
CC-2324.5-NS	23" X 24.5"
CC-2426-NS**	24" X 26"

Nitrogen Tank Covers - 50 p/case

Product	Dimensions
NTC-1862-NS	18" X 62"

Case Cart Covers - 100 p/roll

Product	Dimensions
NB-7248-G-NS	72" X 48"
ONB-7248-G-NS	72" X 48" (opaque)

Please allow 5-7 days for custom non-sterile orders to ship!

Ordering Information

Fax orders may be placed 24 hours a day, 7 days a week.

Orders placed after 4 pm Eastern or on weekends
will be processed the next business day.

When ordering please have the following information:

Shipping and Billing Address

Purchase Order Number

Phone Number • All product information

Name of department ordering

**Our fully trained
staff is eager to
assist you!**



Miranda Standridge
Head of Sales Department
mirandas@preferred-medical.com
Ext. 302



Heather Russell
Sales Representative
heatherl@preferred-medical.com
Ext. 303



Ann Fuller
Customer Service
annf@preferred-medical.com
Ext. 402



Taylor Towe
Customer Service
taylort@preferred-medical.com
Ext. 401



Stacy Mullinax
Office Assistant
stacym@preferred-medical.com
Ext. 204



Stephanie Walker
Admin/Office Manager
stephaniew@preferred-medical.com
Ext. 202

**Office Hours:
8 am - 5 pm EST**

**Call (800) 441-1161
Fax (423) 496-3758**

SAME DAY SHIPPING ON MOST ITEMS!



FedEx



Shipping:

- ❖ All shipments are prepay and add via Regular Ground Service unless otherwise requested
- ❖ Normal delivery time is 3-4 business days
- ❖ A variety of special shipping methods are available upon request

Terms:

- ❖ Please pay by invoice date, Net 30 days
- ❖ Payments are to be remitted to:
Preferred Medical Products LLC
PO Box 100
Ducktown, TN 37326



Return and Exchange Policy:

- ♦ No return will be accepted without a Return Goods Authorization (RGA) number received from Preferred Medical Products
- ♦ Returns will be allowed up until 90 days of invoice
- ♦ Custom items are not normally returnable, however, if another customer purchases the item then it may be returned
- ♦ There will be a 15% restocking fee on any return received back to Preferred Medical Products in a sellable condition. No credit will be issued on non-restockable items.
- ♦ Any return due to the fault of Preferred Medical Products will be replaced or credited in full, including freight charges. In this event, a call tag will be issued.
- ♦ Customer will be liable for return freight and original shipping freight on any return due to low usage, overstock, etc.

The logo for Preferred Medical Products LLC features the word "Preferred" in a large, elegant, dark blue script font. Below it, "Medical Products LLC" is written in a smaller, dark blue serif font. To the left of the text is a red diamond shape. Below the company name, the address "191 Industrial Drive • PO Box 100" and "Ducktown, TN 37326" are listed in a dark blue sans-serif font.

Preferred
Medical Products LLC
191 Industrial Drive • PO Box 100
Ducktown, TN 37326



We Appreciate Your Business!

Order via

 Phone (423) 496-3738
(800) 441-1161

 Fax (423) 496-3758

 Online www.preferred-medical.com

Or Email Your Sales Representative!